

Think. Inside. The. Box.

PURELAB® flex 2*



DI Pack

LC208

Lifetime expectancy: 12 months

What is purified from the water?

Ions

Impact of not changing:

Weakly bound ions are released back into the water reducing resistivity



UV light (185nm)

LC210

Lifetime expectancy: 12-18 months

What is purified from the water?

Purifies TOC, inactivates DNA and RNA

Impact of not changing:

Uncontrolled bacterial contamination and TOC levels are not reduced



Point of Use Filter**

LC145

Lifetime expectancy: 3 months

What is purified from the water?

Particulates and bacteria


Impact of not changing:

Potential recontamination of purified water and reduced flow rates

*A combination of treatments and technologies produce the final water specification

**For DNA, RNA & endotoxin specification at point of use, please use a LC197 Biofilter

For more information please contact your local ELGA water specialist:

 (877)-315 3542

 Elga.orders@veolia.com

 e-shop!

