

Think. Inside. The. Box.

PURELAB® flex 4*



UV light (185nm) LC210

Lifetime expectancy: 12-18 months

What is purified from the water?

Purifies TOC and inactivates DNA and RNA

Impact of not changing:

Uncontrolled bacterial contamination and TOC levels are not reduced

DI Pack LC214

Lifetime expectancy: 6 months

What is purified from the water?

Ions

Impact of not changing:

Weakly bound ions are released back into the water reducing resistivity



Composite Vent Filter LC216

Lifetime expectancy: 12 months

What protection does it provide?

Reduces CO₂ and airborne organics from entering the sealed system

Impact of not changing:

Increase in airborne contaminants and potential recontamination of purified

water

Point of Use Filter** LC145

Lifetime expectancy: 3 months

What is purified from the water?

Particulates and bacteria

Impact of not changing:


Potential recontamination of purified water and reduced flow rates



*A combination of treatments and technologies produce the final water specification

**For DNA, RNA & endotoxin specification at point of use, please use a LC197Biofilter

For more information please contact your local ELGA water specialist:

 (877)-315 3542

 Elga.orders@veolia.com

 e-shop!

